

Anti-Dcaf7 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P61962

Item Number: CM0011292

Introduction

High-specificity anti-DCAF7 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat DCAF7. Ideal for craniofacial development and ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	DCAF7 Rabbit mAb
Remarks/Alias	AN11; HAN11; WDR68; SWAN-1; DCAF7
Species	Human
GeneID (Human)	10238
GeneID	10238
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human DCAF7 (P61962).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61962
Protein Weight	39kDa
Shipping	Ice bag