

Anti-Dcaf7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61962

Item Number: CM0017790

Introduction

Rabbit polyclonal antibody against DCAF7, also known as WDR68, for use in Western blot, immunofluorescence, and ELISA. Validated in human and mouse samples. 39kDa protein.

[Learn More](#)

Parameter	Specification
Product Name	DCAF7 Rabbit pAb
Remarks/Alias	AN11; HAN11; WDR68; SWAN-1; DCAF7
Species	Human
GeneID (Human)	10238
GeneID	10238
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-342 of human DCAF7 (NP_005819.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61962
Protein Weight	39kDa
Shipping	Ice bag