

Anti-Dcamkl1/Dclk1 Polyclonal Antibody For Wb, If-P, Elisa - O15075

Item Number: CM0029405



Introduction

Rabbit polyclonal antibody specific for DCAMKL1/DCLK1 (amino acids 590-740). Validated for WB, IF-P, and ELISA in human, mouse, and rat samples. Research-use only for studying neuronal migration and kinase signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	DCAMKL1/DCLK1 Rabbit pAb
Remarks/Alias	CL1; DCLK; CLICK1; DCDC3A; DCAMKL1; DCAMKL1/DCLK1
Species	Human
GeneID (Human)	9201
GeneID	9201
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 590-740 of DCAMKL1/DCLK1 (NP_001317000.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15075
Protein Weight	82kDa
Shipping	Ice bag