

Anti-Dcamkl1/Dclk1 Rabbit Monoclonal Antibody For Wb, Elisa - O15075

Item Number: CM0011902



Introduction

Rabbit monoclonal antibody targeting human DCAMKL1/DCLK1 (serine/threonine-protein kinase), validated for Western blotting and ELISA, with cross-reactivity to mouse and rat. Ideal for neuroscience research on neuronal migration and brain development.

[Learn More](#)

Parameter	Value
Product Name	DCAMKL1/DCLK1 Rabbit mAb
Remarks/Alias	CL1; DCLK; CLICK1; DCDC3A; DCAMKL1; DCAMKL1/DCLK1
Species	Human
GeneID	9201 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 600-700 of human DCAMKL1/DCLK1 (O15075)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15075
Protein Weight	82kDa
Shipping	Ice bag