


Anti-Dck Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P27707

Item Number: CM0007635



Introduction

Accurate rabbit monoclonal antibody against human DCK (deoxycytidine kinase), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, rat. Essential for nucleoside analog phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] DCK Rabbit mAb
Remarks/Alias	DCK; [KD Validated] DCK
Species	Human
GeneID (Human)	1633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DCK (NP_000779.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27707
Protein Weight	31kDa
Shipping	Ice bag