


Anti-Dck Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27707

Item Number: CM0023671



Introduction

Rabbit polyclonal antibody targeting human deoxycytidine kinase (DCK) for WB, IHC-P, IF/ICC, and ELISA. Validated in human, mouse, and rat samples; ideal for cancer and antiviral drug activation studies.

[Learn More](#)

Parameter	Value
Product Name	DCK Rabbit pAb
Remarks/Alias	DCK
Species	Human
GeneID (Human)	1633
GeneID	1633
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DCK (NP_000779.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27707
Protein Weight	31kDa
Shipping	Ice bag