


Anti-Dclre1C Rabbit Polyclonal Antibody For Wb, Elisa - Q96Sd1

Item Number: CM0026487



Introduction

Rabbit polyclonal antibody against human DCLRE1C (Artemis) for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying DNA repair, V(D)J recombination, and immunodeficiency diseases.

[Learn More](#)

Parameter	Value
Product Name	DCLRE1C Rabbit pAb
Remarks/Alias	SCIDA, SNM1C, A-SCID, RS-SCID, DCLREC1C, DCLRE1C
Species	Human
Gene ID	64421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human DCLRE1C (NP_001029027.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96SD1
Protein Weight	78kDa
Shipping	Ice bag