

Anti-Dcn (Decorin) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07585

Item Number: CM0018111

Introduction

High-quality rabbit polyclonal antibody against human DCN/Decorin (P07585), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for extracellular matrix and fibrillogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	DCN Rabbit pAb
Remarks/Alias	CSCD; PG40; PGII; PGS2; DSPG2; SLRR1B; DCN
Species	Human
Gene ID (Human)	1634
Gene ID	1634
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 31-359 of human DCN (NP_001911.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07585
Protein Weight	40kDa
Shipping	Ice bag