


Anti-Dcp1A Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Npi6

Item Number: CM0011237



Introduction

High-quality Anti-DCP1A Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human and mouse. Suitable for mRNA decapping and nonsense-mediated mRNA decay research.

[Learn More](#)

Parameter	Value
Product Name	DCP1A Rabbit mAb
Remarks/Alias	SMIF, SMAD4IP1, HSA275986, Nbla00360, DCP1A
Species	Human
GeneID (Human)	55802
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 483-582 of human DCP1A (Q9NPI6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NPI6
Protein Weight	63kDa
Shipping	Ice bag