

Anti-Dcp1B Polyclonal Antibody For Wb, Elisa - Q8Izd4

Item Number: CM0019266

Introduction

Rabbit polyclonal antibody against human DCP1B, validated for Western blot and ELISA. Cross-reacts with human and rat. Ideal for studying mRNA decapping mechanisms and nonsense-mediated decay.

[Learn More](#)

Parameter	Value
Product Name	DCP1B Rabbit pAb
Remarks/Alias	DCP1; DCP1B
Species	Human
GeneID (Human)	196513
GeneID	196513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 386-617 of human DCP1B (NP_689853.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IZD4
Protein Weight	68kDa
Shipping	Ice bag