

Anti-Dcp2 Rabbit Polyclonal Antibody For Wb, Elisa - Q8lu60

Item Number: CM0029779



Introduction

Rabbit polyclonal antibody targeting DCP2 (NUDT20), validated for WB and ELISA. Reacts with human, mouse, rat. DCP2 is an mRNA decapping enzyme involved in mRNA turnover and nonsense-mediated decay. Ideal for studying post-transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	DCP2 Rabbit pAb
Remarks/Alias	NUDT20; DCP2
Species	Human
Gene ID	167227
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human DCP2 (NP_001229306.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IU60
Protein Weight	48kDa
Shipping	Ice bag