

## Anti-Dcps Polyclonal Antibody For Wb, Elisa - Q96C86

Item Number: CM0013629



### Introduction

Rabbit polyclonal antibody targeting human DCPS (m7GpppX diphosphatase). Suitable for Western blot and ELISA. Cross-reacts with mouse. Detects ~39 kDa band. Essential for mRNA decay and splicing research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | DCPS Rabbit Polyclonal Antibody                                                                                  |
| Remarks/Alias    | ARS; DCS1; HSL1; HINT5; HINT-5; HSPC015; DCPS                                                                    |
| Species          | Human                                                                                                            |
| GeneID           | 28960                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human DCPS (NP_054745.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse                                                                                                     |
| SWISS            | Q96C86                                                                                                           |
| Protein Weight   | 39kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |