

Anti-Dcr1/Cd263 Rabbit Monoclonal Antibody For Wb, Elisa - O14798

Item Number: CM0007522

Introduction

High-specificity rabbit monoclonal antibody against human DcR1/CD263 (TRAIL-R3, TNFRSF10C), validated for Western blot and ELISA. This decoy receptor competes for TRAIL binding to inhibit apoptosis. Reacts with human, 27kDa.

[Learn More](#)

Parameter	Value
Product Name	DcR1/CD263 Rabbit mAb
Remarks/Alias	LIT; DCR1; TRID; CD263; TRAILR3; TRAIL-R3; DCR1-TNFR; DcR1 / CD263
Species	Human
GeneID (Human)	8794
GeneID	8794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-230 of human DcR1 / CD263 (NP_003832.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14798
Protein Weight	27kDa
Shipping	Ice bag