

Anti-Dcr1/Cd263 Rabbit Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - O14798

Item Number: CM0005773

Introduction

High-quality rabbit monoclonal antibody against human DcR1/CD263 (TNFRSF10C, TRAIL-R3) for WB, IHC-P, FC, and ELISA applications. Ideal for apoptosis and TRAIL signaling research. Reliable raw material for IVD development.

[Learn More](#)

Parameter	Value
Product Name	DcR1/CD263 Rabbit mAb
Remarks/Alias	TNFRSF10C; CD263; DCR1; DCR1-TNFR; LIT; TRAIL-R3; TRAILR3; TRID; TNF receptor superfamily member 10c; DcR1/CD263
Species	Human
GeneID (Human)	8794
GeneID	8794
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-235 of human DcR1/CD263 (NP_003832.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	O14798
Protein Weight	27kDa
Shipping	Ice bag