

Anti-Dcr2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubn6

Item Number: CM0009570

Introduction

High-quality anti-DCR2 rabbit monoclonal antibody validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for studying TRAIL receptor signalling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	DCR2 Rabbit mAb
Remarks/Alias	DCR2; CD264; TRUNDD; TRAILR4; TRAIL-R4
Species	Human
Gene ID	8793
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 287-386 of human DCR2 (Q9UBN6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBN6
Protein Weight	42kDa
Shipping	Ice bag