

Anti-Dct Polyclonal Antibody For Wb, Elisa - P40126

Item Number: CM0017108



Introduction

Anti-DCT rabbit polyclonal antibody for WB and ELISA applications. Specific to human and mouse DCT (TRP-2), a 59kDa melanogenic enzyme involved in melanin biosynthesis. Ideal for pigmentation research.

[Learn More](#)

Parameter	Value
Product Name	DCT Rabbit pAb
Remarks/Alias	OCA8; TRP-2; TYRP2; DCT
Species	Human
GeneID (Human)	1638
GeneID	1638
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-350 of human DCT (NP_001913.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40126
Protein Weight	59kDa
Shipping	Ice bag