


Anti-Dctn2 Polyclonal Antibody For Wb, Elisa - Q13561

Item Number: CM0014278



Introduction

Rabbit polyclonal antibody targeting human DCTN2 (dynactin subunit 2). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Suitable for studying dynein-mediated transport and mitosis.

[Learn More](#)

Parameter	Value
Product Name	DCTN2 Rabbit pAb
Remarks/Alias	RBP50; DCTN50; HEL-S-77; DYNAMITIN; DCTN2
Species	Human
GeneID (Human)	10540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-401 of human DCTN2 (NP_006391.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13561
Protein Weight	44kDa
Shipping	Ice bag