

## Anti-Dctn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13561

Item Number: CM0019218



### Introduction

Detect dynactin subunit 2 (DCTN2/p50 dynamitin, 44kDa) in human and rat with this rabbit polyclonal antibody. Validated for Western blot, immunofluorescence/ICC, and ELISA. Ideal for studies on dynein-mediated transport and mitosis.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | DCTN2 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | RBP50; DCTN50; HEL-S-77; DYNAMITIN; DCTN2                                                                                               |
| Species          | Human                                                                                                                                   |
| GeneID (Human)   | 10540                                                                                                                                   |
| GeneID           | 10540                                                                                                                                   |
| Immunogen        | Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 50-320 of human DCTN2 (NP_006391.1)) |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                                       |
| Cross Reactivity | Human, Rat                                                                                                                              |
| SWISS            | Q13561                                                                                                                                  |
| Protein Weight   | 44kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |