


Anti-Dctn5 Monoclonal Antibody For Wb, Elisa - Q9Bte1

Item Number: CM0003341



Introduction

DCTN5 Rabbit mAb, a monoclonal antibody targeting dynactin subunit 5 (p25), validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for dynein transport and kinetochore research.

[Learn More](#)

Parameter	Value
Product Name	DCTN5 Rabbit mAb
Remarks/Alias	DCTN5; dynactin subunit 5
Species	Human
Gene ID (Human)	84516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human DCTN5 (Q9BTE1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTE1
Protein Weight	20kDa
Shipping	Ice bag