


Anti-Dctn5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bte1

Item Number: CM0020909



Introduction

High-quality anti-DCTN5 rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat DCTN5. Ideal for studying dynactin complex function and dynein-mediated intracellular transport.

[Learn More](#)

Parameter	Value
Product Name	DCTN5 Rabbit pAb
Remarks/Alias	DCTN5
Species	Human
GeneID (Human)	84516
GeneID	84516
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DCTN5 (NP_115875.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BTE1
Protein Weight	20kDa
Shipping	Ice bag