

Anti-Dcxr Polyclonal Antibody For Wb And Elisa - Q7Z4W1

Item Number: CM0025949



Introduction

High-quality rabbit polyclonal antibody against human DCXR (L-xylulose reductase). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for uronate cycle and osmoregulation research.

[Learn More](#)

Parameter	Value
Product Name	DCXR Rabbit pAb
Remarks/Alias	XR; DCR; HCR2; P34H; HCR1; KIDCR; PNTSU; SDR20C1; DCXR
Species	Human
Gene ID (Human)	51181
Gene ID	51181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-244 of human DCXR (NP_057370.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z4W1
Protein Weight	26kDa
Shipping	Ice bag