


Anti-Ddah1 Rabbit Polyclonal Antibody For Wb, If/Icc And Elisa - O94760

Item Number: CM0015862



Introduction

High-quality rabbit polyclonal antibody against DDAH1, validated for WB, IF/ICC, and ELISA applications. Ideal for studying nitric oxide regulation and dimethylarginine metabolism in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	DDAH1 Rabbit pAb
Remarks/Alias	DDAH; DDAH1; DDAH-1; HEL-S-16; DDAH1
Species	Human
Gene ID	23576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human DDAH1 (NP_036269.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94760
Protein Weight	31kDa
Shipping	Ice bag