

Anti-Ddah2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95865

Item Number: CM0025767

Introduction

Rabbit polyclonal antibody targeting human DDAH2, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Supports research into DDAH2's roles in vascular biology, insulin secretion, and mitochondrial fission.

[Learn More](#)

Parameter	Value
Product Name	DDAH2 Rabbit pAb
Remarks/Alias	G6a; DDAH; NG30; DDAHII; HEL-S-277; DDAH2
Species	Human
Gene ID (Human)	23564
Gene ID	23564
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human DDAH2 (NP_039268.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95865
Protein Weight	30kDa
Shipping	Ice bag