

Anti-Ddb1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16531

Item Number: CM0027317



Introduction

Anti-DDB1 Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat DDB1 (127kDa). Used for DNA repair and ubiquitination research via UV-DDB and DCX complexes.

[Learn More](#)

Parameter	Value
Product Name	DDB1 Rabbit pAb
Remarks/Alias	XPE; DDBA; XAP1; XPCE; XPE-BF; UV-DDB1; WHIKERS; DDB1
Species	Human
Gene ID (Human)	1642
Gene ID	1642
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 741-1140 of human DDB1 (NP_001914.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16531
Protein Weight	127kDa
Shipping	Ice bag