


Anti-Ddb2 Polyclonal Antibody For If/Icc, Elisa - Q92466

Item Number: CM0019031



Introduction

DDB2 rabbit polyclonal antibody validated for IF/ICC and ELISA in human samples. Targets DNA damage-binding protein 2 (UV-DDB2), a key factor in nucleotide excision repair. Ideal for DNA damage and repair research.

[Learn More](#)

Parameter	Value
Product Name	DDB2 Rabbit pAb
Remarks/Alias	XPE; DDBB; UV-DDB2; DDB2
Species	Human
GeneID	1643
Immunogen	Recombinant fusion protein corresponding to amino acids 1-315 of human DDB2 (NP_000098.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q92466
Protein Weight	48kDa
Shipping	Ice bag