


## Anti-Ddb2 Rabbit Monoclonal Antibody For Wb, Elisa - Q92466

Item Number: CM0003355



### Introduction

This rabbit monoclonal antibody against human DDB2 (DNA damage-binding protein 2, SWISS: Q92466) is validated for Western blot and ELISA. It recognizes a key nucleotide excision repair factor and E3 ligase adaptor, supporting research in genomic integrity and cancer.

[Learn More](#)

| Parameter        | Value                    |
|------------------|--------------------------|
| Product Name     | DDB2 Rabbit mAb          |
| Remarks/Alias    | XPE; DDBB; UV-DDB2; DDB2 |
| Species          | Human                    |
| GeneID (human)   | 1643                     |
| GeneID           | 1643                     |
| Immunogen        | Recombinant protein      |
| Source           | Rabbit                   |
| Category         | Monoclonal Antibodies    |
| Application      | WB, ELISA                |
| Cross Reactivity | Human                    |
| SWISS            | Q92466                   |
| Protein Weight   | 48kDa                    |
| Shipping         | Ice bag                  |