


Anti-Ddb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92466

Item Number: CM0024011



Introduction

Rabbit polyclonal antibody against human DDB2 (DNA damage-binding protein 2), validated for WB, IF/ICC, and ELISA. Suitable for DNA repair research in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	DDB2 Rabbit pAb
Remarks/Alias	XPE; DDBB; UV-DDB2; DDB2
Species	Human
GeneID (Human)	1643
GeneID	1643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-315 of human DDB2 (NP_000098.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92466
Protein Weight	48kDa
Shipping	Ice bag