


Anti-Ddc Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P20711

Item Number: CM0028412



Introduction

High-quality rabbit polyclonal antibody against human DDC, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for catecholamine and serotonin biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	DDC Rabbit pAb
Remarks/Alias	AADC; DDC
Species	Human
GeneID (human)	1644
GeneID	1644
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-137 of human DDC (NP_000781.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20711
Protein Weight	54kDa
Shipping	Ice bag