

Anti-Ddit3/Chop Monoclonal Antibody For Wb, Ihc-P, Elisa - P35638

Item Number: CM0003720

Introduction

Rabbit monoclonal antibody targeting human DDIT3/CHOP (DNA damage-inducible transcript 3). Validated for WB, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Predicted molecular weight 19kDa.

[Learn More](#)

Parameter	Value
Product Name	DDIT3/CHOP Rabbit mAb
Remarks/Alias	CHOP; CEBPZ; CHOP10; CHOP-10; GADD153; AltDDIT3; C/EBPzeta; DDIT3/CHOP
Species	Human
GeneID (Human)	1649
GeneID	1649
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human DDIT3/CHOP (NP_004074.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35638
Protein Weight	19kDa
Shipping	Ice bag