

Anti-Ddit3/Chop Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35638

Item Number: CM0004576

Introduction

High-quality rabbit monoclonal antibody against human DDIT3/CHOP, validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat. Ideal for ER stress and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DDIT3/CHOP Rabbit mAb
Remarks/Alias	CHOP; CEBPZ; CHOP10; CHOP-10; GADD153; AltDDIT3; C/EBPzeta; DDIT3/CHOP
Species	Human
GeneID (Human)	1649
GeneID	1649
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DDIT3/CHOP (NP_004074.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35638
Protein Weight	19 kDa
Shipping	Ice bag