

Anti-Ddit3/Chop Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P35638

Item Number: CM0015894

Introduction

Rabbit polyclonal antibody targeting DDIT3/CHOP (CHOP10, GADD153). Validated in WB, IHC-P, IP, and ELISA across human, mouse, and rat samples. Ideal for investigating ER stress, apoptosis, and unfolded protein response pathways.

[Learn More](#)

Parameter	Value
Product Name	DDIT3/CHOP Rabbit pAb
Remarks/Alias	CHOP; CEBPZ; CHOP10; CHOP-10; GADD153; AltDDIT3; C/EBPzeta; DDIT3/CHOP
Species	Human
Gene ID (Human)	1649
Gene ID	1649
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human DDIT3/CHOP (NP_004074.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35638
Protein Weight	19kDa
Shipping	Ice bag