

Anti-Ddost Polyclonal Antibody For Wb, Ihc-P, Elisa - P39656

Item Number: CM0029500



Introduction

DDOST rabbit polyclonal antibody suitable for WB, IHC-P, ELISA. Detects human, mouse, and rat DDOST, a key OST subunit essential for protein N-glycosylation. Ideal for glycobiology research.

[Learn More](#)

Parameter	Value
Product Name	DDOST Rabbit pAb
Remarks/Alias	OST; WBP1; AGER1; CDG1R; GATD6; OST48; OKSWcl45; DDOST
Species	Human
Gene ID (Human)	1650
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human DDOST (NP_005207.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P39656
Protein Weight	51kDa
Shipping	Ice bag