

Anti-Ddr1/Mck10/Cd167A Rabbit Monoclonal Antibody For Wb And Elisa - Q08345

Item Number: CM0006780



Introduction

Rabbit monoclonal antibody against human DDR1/MCK10/CD167a, validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying collagen receptor signaling, cell migration, and tumor invasion pathways.

[Learn More](#)

Parameter	Value
Product Name	DDR1/MCK10/CD167a Rabbit mAb
Remarks/Alias	CAK; DDR; NEP; HGK2; PTK3; RTK6; TRKE; CD167; EDDR1; MCK10; NTRK4; PTK3A
Species	Human
GeneID (human)	780
GeneID	780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-416 of human DDR1/MCK10/CD167a (NP_054699.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q08345
Protein Weight	101kDa
Shipping	Ice bag