


Anti-Ddr2 Rabbit Monoclonal Antibody For Wb And Elisa - Q16832

Item Number: CM0005378



Introduction

High-quality rabbit monoclonal anti-DDR2 antibody validated for Western blot and ELISA applications. Detects human, mouse, and rat DDR2, a collagen receptor tyrosine kinase involved in tissue remodeling and bone development.

[Learn More](#)

Parameter	Value
Product Name	DDR2 Rabbit mAb
Remarks/Alias	TKT; WRCN; MIG20a; NTRKR3; TYRO10; DDR2
Species	Human
GeneID (Human)	4921
GeneID	4921
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 756-855 of human DDR2 (Q16832).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16832
Protein Weight	97kDa
Shipping	Ice bag