


Anti-Ddr2 Rabbit Polyclonal Antibody For Wb, Elisa - Q16832

Item Number: CM0012945



Introduction

Rabbit polyclonal antibody against human DDR2, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Target involved in collagen signaling and tissue remodeling; suitable for oncology and bone development research.

[Learn More](#)

Parameter	Specification
Product Name	DDR2 Rabbit pAb
Remarks/Alias	TKT; WRCN; MIG20a; NTRKR3; TYRO10; DDR2
Species	Human
Gene ID (Human)	4921
Gene ID	4921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-400 of human DDR2 (NP_001014796.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16832
Protein Weight	97kDa
Shipping	Ice bag