


Anti-Ddt Rabbit Polyclonal Antibody For Wb, Elisa - P30046

Item Number: CM0027223



Introduction

Rabbit polyclonal antibody against human D-dopachrome decarboxylase (DDT). Validated for Western blotting and ELISA with cross-reactivity to human, mouse, and rat. Ideal for melanin biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	DDT Rabbit pAb
Remarks/Alias	D-DT; DDCT; MIF2; MIF-2; DDT
Species	Human
GeneID (Human)	1652
GeneID	1652
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human DDT (NP_001346.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30046
Protein Weight	13kDa
Shipping	Ice bag