

Anti-Ddx1 Polyclonal Antibody - Q92499

Item Number: CM0027325



Introduction

Rabbit polyclonal DDX1 antibody for WB, IHC-P, IF/ICC, IP, ELISA. Validated in human, mouse, rat. Detects ATP-dependent RNA helicase DDX1 involved in transcription, RNA processing, and antiviral defense - Q92499.

[Learn More](#)

Parameter	Value
Product Name	DDX1 Rabbit pAb
Remarks/Alias	DBP-RB; UKVH5d; DDX1
Species	Human
GeneID (Human)	1653
GeneID	1653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 441-740 of human DDX1 (NP_004930.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92499
Protein Weight	82kDa
Shipping	Ice bag