

Anti-Ddx17 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92841

Item Number: CM0011944



Introduction

DDX17 rabbit monoclonal antibody validated for WB, IHC-P, ELISA with cross-reactivity in human, mouse, rat. Targets RNA helicase DDX17 (Q92841) involved in pre-mRNA splicing and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	DDX17 Rabbit mAb
Remarks/Alias	P72; RH70; DDX17
Species	Human
GeneID (Human)	10521
GeneID	10521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human DDX17 (Q92841)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92841
Protein Weight	80kDa
Shipping	Ice bag