


Anti-Ddx21 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nr30

Item Number: CM0016286



Introduction

Rabbit polyclonal antibody targeting human DDX21 (Nucleolar RNA helicase 2). Validated for WB, IHC-P, IF/ICC, ELISA. Detects 87 kDa protein. Cross-reacts with human and mouse. Ideal for rRNA processing, transcription regulation, and antiviral innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	DDX21 Rabbit pAb
Remarks/Alias	RH; GUA; GURDB; II/Gu; RH II/Gu; RH-II/GU; gu-alpha; RH-II/GuA; DDX21
Species	Human
Gene ID (Human)	9188
Gene ID	9188
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NR30
Protein Weight	87kDa
Shipping	Ice bag