

Anti-Ddx24 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9GZr7

Item Number: CM0018543

Introduction

Rabbit polyclonal antibody to human DDX24 (DEAD box protein 24), an ATP-dependent RNA helicase involved in pre-mRNA splicing and innate immunity. Validated for WB, IF/ICC, ELISA in human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	DDX24 Rabbit pAb
Remarks/Alias	DDX24
Species	Human
GeneID (human)	57062
GeneID	57062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-859 of human DDX24 (NP_065147.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZR7
Protein Weight	96kDa
Shipping	Ice bag