

Anti-Ddx3 Rabbit Polymab® Antibody For Wb, If/Icc, Ip, Elisa - O00571

Item Number: CM0006865



Introduction

High-quality rabbit monoclonal antibody against DDX3 (ATP-dependent RNA helicase DDX3X) validated in WB, IF/ICC, IP, ELISA for human, mouse, rat. Ideal for transcription, translation, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	DDX3 Rabbit PolymAb®
Remarks/Alias	DBX; DDX3; HLP2; DDX14; CAP-Rf; MRX102; MRXSSB
Species	Human
GeneID	1654
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human DDX3 (NP_001347.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00571
Protein Weight	73kDa
Shipping	Ice bag