

Anti-Ddx31 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H8H2

Item Number: CM0024360

Introduction

DDX31 Rabbit pAb is a polyclonal antibody for WB, IF/ICC, ELISA, reacting with human, mouse, rat. Targets the DEAD box helicase DDX31, essential for ribosome biogenesis and p53 regulation.

[Learn More](#)

Parameter	Value
Product Name	DDX31 Rabbit pAb
Remarks/Alias	PPP1R25; DDX31
Species	Human
Gene ID (Human)	64794
Gene ID	64794
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 652-851 of human DDX31 (NP_073616.6).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8H2
Protein Weight	94kDa
Shipping	Ice bag