


Anti-Ddx39A Rabbit Polyclonal Antibody - 000148

Item Number: CM0030473



Introduction

Rabbit polyclonal antibody targeting human DDX39A (BAT1/URH49), validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for studying mRNA export, splicing, and antiviral mechanisms.

[Learn More](#)

Parameter	Value
Product Name	DDX39A Rabbit pAb
Remarks/Alias	BAT1; DDXL; BAT1L; DDX39; URH49; DDX39A
Species	Human
GeneID (Human)	10212
Immunogen	Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DDX39A (NP_005795.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00148
Protein Weight	49kDa
Shipping	Ice bag