


Anti-Ddx3X Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O00571

Item Number: CM0024798



Introduction

Rabbit polyclonal antibody against human DDX3X (O00571). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with Mouse, Rat. ATP-dependent RNA helicase involved in transcription, translation, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	DDX3X Rabbit pAb
Remarks/Alias	DBX; DDX3; HLP2; DDX14; CAP-Rf; MRX102; MRX55B; DDX3X
Species	Human
GeneID (Human)	1654
GeneID	1654
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DDX3X (NP_001347.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00571
Protein Weight	73kDa
Shipping	Ice bag