


Anti-Ddx3X/Ddx3Y Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O00571/O15523

Item Number: CM0007009



Introduction

High-quality rabbit monoclonal antibody against DDX3X/DDX3Y, validated for Western blot, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying RNA helicase functions, viral infections, and innate immunity signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	DDX3X/DDX3Y Rabbit mAb
Remarks/Alias	DBX; DDX3; HLP2; DDX14; CAP-Rf; MRX102; MRXSSB; DBY; DDX3X/DDX3Y
Species	Human
Gene ID (Human)	1654/8653
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00571/O15523
Protein Weight	73 kDa/71 kDa/32 kDa
Shipping	Ice bag