

Anti-Ddx3Y Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15523

Item Number: CM0017072

Introduction

High-affinity anti-DDX3Y polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat DDX3Y protein. Ideal for studying ATP-dependent RNA helicase function and immune response pathways.

[Learn More](#)

Parameter	Value
Product Name	DDX3Y Rabbit pAb
Remarks/Alias	DBY; DDX3Y
Species	Human
GeneID (human)	8653
GeneID	8653
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human DDX3Y (NP_001116137.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15523
Protein Weight	73kDa
Shipping	Ice bag