

Anti-Ddx3Y Rabbit Polyclonal Antibody For Wb, Elisa - 015523

Item Number: CM0016877

Introduction

Rabbit polyclonal antibody targeting human DDX3Y, validated for WB and ELISA. Reacts with human and rat. Suitable for studying RNA helicase function in immune response and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	DDX3Y Rabbit pAb
Remarks/Alias	DBY; DDX3Y
Species	Human
Gene ID (Human)	8653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human DDX3Y (NP_001116137.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O15523
Protein Weight	73kDa
Shipping	Ice bag