

Anti-Ddx4 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Nqi0

Item Number: CM0014562

Introduction

Rabbit polyclonal antibody against DDX4 (VASA), a germ cell-specific ATP-dependent RNA helicase. Detects human, mouse, and rat DDX4 for Western blot, immunofluorescence, IHC, and ELISA applications in spermatogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DDX4 Rabbit pAb
Remarks/Alias	VASA; DDX4
Species	Human
Gene ID	54514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 635-724 of human DDX4 (NP_077726.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NQI0
Protein Weight	79kDa
Shipping	Ice bag