

Anti-Ddx4 (Vasa) Rabbit Monoclonal Antibody For Wb - Q9Nqi0

Item Number: CM0004162

Introduction

High-quality anti-DDX4 (VASA) rabbit monoclonal antibody validated for WB, IP, IF-P, IHC-P, and ELISA; cross-reacts with human, mouse, rat. Ideal for germ cell development and spermatogenesis research.

[Learn More](#)

Parameter	Value
Product Name	DDX4 Rabbit mAb
Remarks/Alias	VASA
Species	Human
Gene ID (Human)	54514
Gene ID	54514
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 635-724 of human DDX4 (NP_077726.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQI0
Protein Weight	79kDa
Shipping	Ice bag