


Anti-Ddx41 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9Ujv9

Item Number: CM0030620



Introduction

Rabbit polyclonal antibody against human DDX41, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Detects human, mouse, and rat DDX41 (70 kDa). Suitable for innate immunity and RNA metabolism research.

[Learn More](#)

Parameter	Value
Product Name	DDX41 Rabbit pAb
Remarks/Alias	ABS; MPLPF; DDX41
Species	Human
GeneID (Human)	51428
GeneID	51428
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human DDX41 (NP_057306.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJV9
Protein Weight	70kDa
Shipping	Ice bag